


Evaluation of Different Temperatures in Cold Air Cooling With Pulsed-Dye Laser

Stefan Hammes and Christian Raulin, MD* Laserklinik Karlsruhe, Karlsruhe, Germany

Table 1
Treatment Results
Mean Values
Scales 0-3
(0 = no, 1 = Little,
2 = Medium, 3 = High)



Background and Objectives

Cold air cooling is widely used in dermatological laser therapy. We investigated the influence of cold air cooling at different skin temperatures on therapeutic outcome and side effects of pulsed dye laser treatment of facial telangiectasia.

Study Design / Materials and Methods

From September 2002 to February 2003, 17 patients with previously untreated facial telangiectasia underwent a single treatment session with flash-lamp pulsed dye laser (3.5 J / cm², 585 nm, 0.45 milliseconds pulse length, 10 mm beam diameter, Cynosure 1 V). The treatment area was divided into three sub-areas: no cooling, cold air cooling to 20°C and to 17°C skin temperature. The skin temperature was monitored by a prototype infrared sensor system which controlled the temperature of the cold air stream (Cryo 5). In a prospective study, we collected data on purpura, pain, clearance, and patient satisfaction on numerical analog scales (NAS) from 0 (meaning »no«) to 3 (meaning »high«).

Results

Without cooling, purpura (2.53), pain (2.41), and clearance (2.35) were rated medium to high. Cooling to 20°C reduced purpura (1.12) and pain (1.06), whereas the clearance (2.12) was only slightly affected. Cooling to 17°C reduced purpura (0.88) and pain (0.76) even more, the clearance (2.06) was lowered marginally. Most patients preferred cooling to 20°C skin temperature.

Conclusion

In dermatological laser therapy of facial telangiectasia, the use of cold air cooling can significantly reduce side effects and increase patient satisfaction while only slightly affecting clearance. Cooling to 20°C skin temperature proved to be a well-balanced middle course. For the practical use of cold air cooling, we thus recommend cooling to a level which the patient can tolerate without problems and to try to increase the energy densities.

Cryo

Technical Data

Power supply	230 V / 50 Hz 230 V / 60 Hz 115 V / 60 Hz
Power input max.	1 KW
Stand-by function	260 W / h
Protection according to IEC 601-1	Class I, Type B
MDD / MPG	Class IIa
Treatment tube length	180 cm
Housing dimensions	H 645 mm / W 390 mm / L 680 mm
Weight	60 kg
Therapy air flow	9 levels, max. 1000 l / min

6 programs combining air flow and treatment time
3 user defined programs
1 favorite user defined program

Set up Menu
Service Menu

Options



Articulating arm



Standard treatment hose



Wheel with brake

Zimmer Elektromedizin GmbH
Junkersstraße 9
D-89231 Neu-Ulm
Tel. +49. (0)7 31. 97 61-0
Fax +49. (0)7 31. 97 61-118
info@zimmer.de
www.zimmer.de

Zimmer
MedizinSysteme

Zimmer
MedizinSysteme



Cryo

Cryo

-30°C Cold Air



**The skin cooling system
d esigned for superficial
laser skin procedures.**
The Cryo 6 Cold Air Device
is intended to minimize
pain and thermal injury
during laser and dermato-
logical treatments and for
temporary topical anesthe-
tic relief for injections.

Unlike other cooling
methods, such as contact
cooling, cryogen spray or
ice packs, the Cryo 6 can
cool the epidermis before,
during and after the laser
energy has been applied,
without deflecting the laser
beam.



Easy to operate
A large display clearly
indicates all treatment
parameters. Cryo 6 fea-
tures a sensory glass
keyboard, which enables
the selection from 6
preset programs.

Select a program. Press
Start. Enjoy the pleasures
waiting at your fingertips.

An optional articulating
arm permits »hand free«
operation.

The light application hose
can be connected to many
laser handpieces. The user
can easily regulate the
air flow for his personal
needs.



User Defined Programs
3 Supplementary storage
possibilities are free for
individual programs.

Individual, Favorite
and Sequence Program
Storage – 3 User-oriented
features providing utmost
treatment efficiency and
state-of-the-art precision.

Storage possibility is
available for the user's
favorite program. This
program appears when
Cryo 6 is turned on and
at the conclusion of each
treatment.

Additionally, the user can
choose 2 sequences to be
stored as a program for
easy selection.

Increased Patient Comfort
Hand-free Operation
No Consumables



Economic
Room temperature air is
filtered and cooled down
to -30°C by a closed
cooling circuitry.
- Cost efficient: no con-
sumable handling or
additional costs
- Powerful: full day oper-
ation with no downtime
- Practical: a special
designed glass shelf just
where you need it – for
a laser, smoke evacuator
or accessories.

Easy Maintenance
A constant monitoring
system for the defrosted
water level and an auto-
matic defrosting feature
provide a smooth-running
daily operation. The air
filter is easy to access, just
vacuum when dusty.