

INTRAcEL PRO before and after

Acne Scars



F 27Y, 6 months post 3rd treatment
Courtesy of Wanitphakdeedecha, MD

Striae



F 26Y, 1 month post 2nd treatment
Courtesy of Greg Chernoff, MD

Acne/Acne Scars



M 42Y, 7 days post treatment
Courtesy of Takashi Takahashi, M.D.

Fine Lines & Wrinkles



F 58Y, 5 weeks post treatment
Courtesy of JD McCoy, NMD



F 40Y, 4 weeks post 2 treatments
Courtesy of Jeisys Medical

Technical Specification

Applied Energy	RF (Radio Frequency) - 1Mhz
Total RF Energy	Up to 78.4 J/cm ²
Spot Size	1 cm x 1cm
Disposable Tip Matrix	49 electrode needles
Trix Energy Array	Up to 1600 mJ MTZ

System Specification

Display	Color LCD Touch Screen
Weight	138 lb. / 63Kg.
Dimensions	13.7 in. x 15.9 in. x 42.5 in.
Electrical Requirements	AC 110 V, 50~60 Hz, 1 KV



Exclusive distributor of Jeisys Medical in the U.S.A.

Perigee Medical - 844-PERIGEE - perigee.com - INTC 09122021 - © Perigee Medical, All rights reserved.

INTRAcEL PRO™


World Leader in MicroRF™ Therapy



The Secret of my Skin

MAXIMUM RESULTS


INTRAcel PRO is the only RF microneedle system that treats 100% of the area uniformly for maximum comfort and results.



INTRAcel PRO BENEFITS

Consistent, more predictable results:

- Precisely delivers energy to the target area
- Controllable depth and intensity of treatment
- Active impedance feedback for optimal RF energy delivery
- Greater epidermal protection, reduced risk of PIH
- Minimal downtime for a high degree of patient satisfaction
- Safe on all skin types
- Very comfortable treatments




INTRAcel PRO™
World Leader in MicroRF Therapy

MicroRF Therapy


OPTIMIZED TIP

(based on 6 years of experience)

INTRAcel PRO tips are optimized to treat various depths of the skin. Each tip has 49 pins specifically designed to reach a selectable depth in the dermis.



The microneedles are insulated except 0.3 mm at the tip where the RF energy is emitted for maximum patient comfort and targeted energy delivery.



Only INTRAcel PRO has patented technology* for precise and painless insertion of insulated needles.

Only INTRAcel PRO offers monopolar energy for full field dermal heating.

Only INTRAcel PRO monitors impedance and assures that optimum treatment parameters are delivered.



*Painless insertion of insulated RF microneedle is proprietary technology protected by U.S. Patent 8,666,487.

"INTRAcel PRO delivers a fractional RF treatment which bypasses the dermal-epidermal junction so we easily minimize risk of PIH and use it on all skin types."
- Michael H. Gold, M.D.

*As reported on the Aesthetic Guide